

Primerdesign

R01059

# Equine Strangles

**Kit version: 1**

## Target regions

*S. equi* subsp. *equi* (eqbE\_3 gene)

*S. equi* subsp. *zooepidemicus*  
(fabG\_2 gene)

genesig<sup>®</sup> PLEX Kit

100 tests

GENESIG

Kits by Primerdesign

For general laboratory and research use only

# Product Description

The genesig®PLEX qPCR kit detects the eqbE\_3 gene from *Streptococcus equi* subspecies *equi* and the fabG\_2 gene from *Streptococcus equi* subspecies *zooepidemicus*.

*Streptococcus equi* are gram-positive, coccoid bacterium with a genome of approximately 2.3 Mb. There are two subspecies:

- *Streptococcus equi* subspecies *equi* (also known as *S. equi* subsp. *equi*, or *S. equi equi*)
- *Streptococcus equi* subspecies *zooepidemicus* (also known as *S. equi* subsp. *zooepidemicus*, or *S. zooepidemicus*)

They typically appear in pairs or long chains under a microscope. They are catalase-negative, an important test to differentiate streptococci from staphylococci and represent as beta-haemolytic when cultured on blood agar. Determination can be difficult as colony morphology between the subspecies can be identical. They are characterised as Lancefield group C and considered to rarely cause disease in humans.

*S. equi equi* is an infectious, contagious disease that is highly adapted to and produces clinical disease in the Equidae family (horses, donkeys and mules). It is highly contagious and produces high morbidity and low mortality in susceptible populations. It is characterised by the formation of abscesses in the lymph tissue of the upper respiratory tract. The formation of abscesses and inflammation that occurs restricts the airways, causing difficulty breathing and swallowing. The resulting strained respiratory noises are responsible for the common name 'strangles'. Symptoms of 'strangles' include high fever, nasal discharge, and abscesses on other areas of the body (Bastard Strangles). Transmission occurs via expelled fomites and direct contact with infectious material such as mucus. Clinically ill horses should be isolated to prevent the spread of the pathogen.








*S. zooepidemicus* is closely related to *S. equi equi*, with genome similarity of ~90%. Studies have indicated that *S. equi equi* is derived from *S. zooepidemicus*. It is an opportunistic pathogen that is a normal part of equine bacterial flora. It is responsible for a wide variety of diseases in susceptible horses including, pneumonia, abortion, chronic skin wounds and epididymitis (the structure that collects and stores sperm). Studies have shown that zoocins released by *S. zooepidemicus* can kill *S. equi equi*. Therefore strangles abscesses that rupture can quickly become colonised by *S. zooepidemicus*.

## Specificity

The genesig®PLEX kit for Equine Strangles is designed for the in vitro detection of *Streptococcus equi* subsp. *equi* and *Streptococcus equi* subsp. *zooepidemicus*. Specifically, the primers will detect over 95% of sequences available on the NCBI database at the time of last review.

The dynamics of genetic variation means that new sequence information may become available after the initial design. If you require further information or have a specific question about the detection profile of this kit then please send an e-mail to [techsupport@primerdesign.co.uk](mailto:techsupport@primerdesign.co.uk) and our team will answer your question.

# Kit contents

Quantity	Component	Tube	Cap Colour
2	<b>Equine Strangles primer/probe mix (50 reactions)</b> S.equi equi (FAM) S.equi zooepidemicus (ROX) Internal Extraction Control (VIC)		<b>BROWN</b>
1	<b>Equine Strangles positive control template</b>		<b>RED</b> (in silver foil wrapper)
2	<b>oasig® Lyophilised 2X qPCR Master Mix (50 reactions per glass vial)</b>		<b>SILVER</b>
2	<b>oasig® Master Mix resuspension buffer</b>		<b>BLUE</b>
1	<b>genesig® Easy DNA internal extraction control</b>		<b>BLUE</b> (in gold foil wrapper)
1	<b>Template preparation buffer</b> for resuspension of internal control template and positive control template		<b>YELLOW</b>
1	<b>RNase/DNase free water</b> for resuspension of primer/probe mixes		<b>WHITE</b>

## Reagents and equipment to be supplied by the user

### Real-time PCR Instrument

Must be able to read fluorescence through FAM, ROX, and HEX/VIC.

### Extraction kit

This kit is recommended for use with genesig® Easy DNA/RNA extraction kit or exsig® Mag. However, it is designed to work well with all processes that yield high-quality nucleic acid with minimal PCR inhibitors.

### Pipettors and filter tips

### Vortex and centrifuge

### 1.5 ml microtubes

### qPCR plates or reaction tubes

## Kit storage and stability

This kit is stable for shipping at ambient temperature but should be stored at -20°C upon arrival. Once the lyophilised components have been resuspended, they should not be exposed to temperatures above -20°C for longer than 30 minutes at a time and unnecessary repeated freeze/thawing should be avoided. The kit is stable for six months from the date of resuspension under these circumstances.

Primer Design Ltd does not recommend using the kit after the expiry date stated on the pack.

## Suitable sample material

This kit can be used with all types of samples from various origins. Please ensure that the extracted nucleic acid samples are suitable in terms of purity, concentration, and DNA/RNA integrity. An internal extraction control is provided to test for non-specific PCR inhibitors.

## Dynamic range of test

Under optimal PCR conditions the kit can achieve priming efficiencies between 90-110% and detect less than 100 copies of target template.

# Principles of the test

## Real-time PCR

Individual primers and probes designed for each target have been combined into a single reaction and these can be detected through different fluorescent channels as described in the kit contents.

The primer and probe mix provided exploits the so-called TaqMan® principle. During PCR amplification, forward and reverse primers hybridize to the target DNA. A fluorogenic probe is included in the same reaction mixture which consists of a DNA probe labelled with a 5'-dye and a 3'-quencher. During PCR amplification, the probe is cleaved, and the reporter dye and quencher are separated. The resulting increase in fluorescence can be detected on a range of qPCR platforms.

## Positive control

For a positive control, the kit contains a single positive control that contains templates for the 2 targets in the test. The kit positive control will give an *S. equi equi* signal through the FAM channel and an *S. equi zooepidemicus* signal through the ROX channel. Each time the kit is used, at least one positive control reaction must be included in the run. A positive result indicates that the primers and probes for detecting each target are working properly in that particular run. If a negative result is obtained the test results are invalid and must be repeated. Care should be taken to ensure that the positive control does not contaminate any other kit component which would lead to false positive results. This can be achieved by handling this component in a post PCR environment. Care should also be taken to avoid cross-contamination of other samples when adding the positive control to the run. This can be avoided by sealing all other samples and negative controls before pipetting the positive control into the positive control well.

## Negative control

To validate any positive findings a negative control reaction should be included every time the kit is used. For this reaction the RNase/DNase free water should be used instead of template. A negative result indicates that the reagents have not become contaminated while setting up the run. It is also known as a No Template Control or NTC.

## Internal DNA extraction control

When performing DNA extraction, it is often advantageous to have an exogenous source of DNA template that is spiked into the lysis buffer. This control DNA is then co-purified with the sample DNA and can be detected as a positive control for the extraction process. Successful co-purification and qPCR for the control DNA also indicates that PCR inhibitors are not present at a high concentration. Within the Equine Strangles primer/probe mix are primers and probes to detect the exogenous DNA using qPCR. The primers are present at PCR limiting concentrations which allows multiplexing with the target sequence primers. Amplification of the control DNA does not interfere with detection of the target DNA even when present at low copy number. The Internal control is detected through the VIC channel and gives a Cq value of 28+/-3.

# Resuspension protocol

To minimise the risk of contamination with foreign DNA, we recommend that all pipetting is performed in a PCR clean environment. Ideally this would be a designated PCR lab or PCR cabinet. Filter tips are recommended for all pipetting steps.

**1. Pulse-spin each tube in a centrifuge before opening.**

This will ensure lyophilised primer/probe mix or template is in the base of the tube and is not lost upon opening the tube.

**2. Resuspend the primer/probe mix in the RNase/DNase free water supplied, according to the table below:**

To ensure complete resuspension allow primer/probe mix to rehydrate for 10 minutes at room temperature. Vortex the tube thoroughly, followed by pipetting up and down 10 times. Failure to mix well can produce poor kit performance.

Component - resuspend in water	Volume
<b>Pre-PCR pack</b>	
Equine Strangles primer/probe mix ( <b>BROWN</b> )	55 µl

**3. Resuspend the internal control template and positive control template in the template preparation buffer supplied, according to the table below:**

To ensure complete resuspension, vortex each tube thoroughly.

Component - resuspend in template preparation buffer	Volume
<b>Pre-PCR heat-sealed foil</b>	
genesig® Easy DNA internal extraction control ( <b>BLUE</b> )	600 µl
<b>Post-PCR heat-sealed foil</b>	
Positive control template ( <b>RED</b> ) *	500 µl

\* This component contains high copy number template and is a VERY significant contamination risk. It must be opened and handled in a separate laboratory environment, away from the other components.

**4. Resuspend the lyophilised oasisig® Master Mix in oasisig® resuspension buffer, according to the table below:**

Component - resuspend in oasisig® resuspension buffer	Volume
oasisig® Lyophilised 2X qPCR Master Mix ( <b>SILVER</b> )	525 µl

# DNA extraction

The internal extraction control DNA can be added either to the DNA lysis/extraction buffer or to the sample once it has been resuspended in lysis buffer.

**DO NOT add the internal extraction control DNA directly to the unprocessed biological sample as this will lead to degradation and a loss in signal.**

1. Add 4 µl of the Internal extraction control DNA (**BLUE**) to each sample in DNA lysis/extraction buffer.
2. Complete nucleic acid extraction according to the manufacturer's protocols.

## qPCR detection protocol

**For optimum performance and sensitivity.**

All pipetting steps and experimental plate set up should be performed on ice. After the plate is prepared proceed immediately to the amplification protocol. Prolonged incubation of reaction mixes at room temperature can lead to PCR artifacts that reduce the sensitivity of detection.

1. For each DNA sample prepare a reaction mix according to the table below:  
Include sufficient reactions for all samples, positive and negative controls.

Component	Volume
oasig <sup>®</sup> Lyophilised 2X qPCR Master Mix ( <b>SILVER</b> )	10 µl
Equine Strangles primer/probe mix ( <b>BROWN</b> )	1 µl
RNase/DNase free water ( <b>WHITE</b> )	4µl
<b>Final Volume</b>	<b>15 µl</b>

2. Pipette 15 µl of this mix into each well according to your experimental qPCR plate set-up.
3. Pipette 5 µl of extracted sample into the sample well according to your experimental plate set-up.  
For negative control wells use 5 µl of RNase/DNase free water (**WHITE**). The final volume for each well is 20 µl.
4. Pipette 5µl of positive control template into the positive control well, according to your experimental plate set up.  
The positive control contains templates for *S. equi equi* and *S. equi zooepidemicus*. The final volume in each well is 20 µl.

# qPCR amplification protocol

Amplification conditions for oasisig<sup>®</sup> Lyophilised 2X qPCR Master Mix (**SILVER**)

	Step	Time	Temp
	Enzyme activation	2 min	95 °C
Cycling x50	Denaturation	10 s	95 °C
	DATA COLLECTION *	60 s	60 °C


\* Fluorogenic data should be collected during this step through the FAM, ROX and VIC channels.




## Interpretation of results

Target	Internal control (VIC)	Positive control	Negative control	Interpretation
FAM+	+ / -	+	-	<b>S. equi equi POSITIVE RESULT</b>
ROX+	+ / -	+	-	<b>S. equi zooepidemicus POSITIVE RESULT</b>
-	+	+	-	<b>NEGATIVE RESULT</b>
+ / -	+ / -	+	$\leq 35$	<b>EXPERIMENT FAILED</b> due to test contamination
+ / -	+ / -	+	$> 35$	*
-	-	+	-	<b>SAMPLE PREPARATION FAILED</b>
+ / -	+ / -	-	+ / -	<b>EXPERIMENT FAILED</b>

\* Where the test sample is positive and the negative control is positive with a Cq  $> 35$ , the sample must be reinterpreted based on the relative signal strength of the two results:



If the sample amplifies  $> 5$  Cq earlier than the negative control, then the sample should be reinterpreted (via the table above) with the negative control verified as negative.



If the sample amplifies  $< 5$  Cq earlier than the negative control, then the positive sample result is invalidated, and the result should be determined inconclusive due to test contamination. The test for this sample should be repeated.

## Positive Control

The positive control well should give an amplification plot through the FAM channel (*S. equi equi*) and the ROX channel (*S. equi zooepidemicus*). There is no internal control template within the positive control so the VIC channel should give no signal (flat amplification plot). The positive control signals indicate that the kit is working correctly to detect each target.

## No template control (NTC)

The NTC should give a flat line (flat amplification plots) through all channels. Signals in the NTC indicate cross contamination during plate set up.

## DNA Internal extraction control

The Cq value obtained with the internal control will vary significantly depending on the extraction efficiency, the quantity of nucleic acid in the PCR reaction and the individual machine settings. Cq values of  $28 \pm 3$  are within the normal range. When amplifying a target sample with a high genome copy number, the internal extraction control may not produce an amplification plot. This does not invalidate the test and should be interpreted as a positive experimental result.

## Sample data

Presence of the target pathogens is shown by amplification in the channels indicated in the kit contents section. Positive signals indicate positive tests for those targets. It may be possible for samples to contain multiple targets, therefore positive results in the FAM and ROX channels may be present.

# Notices and disclaimers

This product is developed, designed and sold for research purposes only. It is not intended for human diagnostic or drug purposes or to be administered to humans unless clearly expressed for that purpose by the Food and Drug Administration in the USA or the appropriate regulatory authorities in the country of use. During the warranty period, Primer Design Ltd genesig® detection kits allow precise and reproducible data recovery combined with excellent sensitivity. For data obtained in violation of the general GLP guidelines and the manufacturer's recommendations, the right to claim under guarantee is expired.

## Trademarks

genesig® is a registered trademark of Primer Design Ltd.

oasig® is a trademark of Primer Design Ltd.

exsig®Mag is a trademark of Primer Design Ltd.