

Primerdesign™ Ltd

# JAK2 V617F mutation

mutation detection by quantitative allele specific  
amplification (quasa)

JAK2(V617F)

50 tests

GENESIG

Kits by Primerdesign

For general laboratory and research use only

# Contents

|  |           |
|--|-----------|
| <b>Introduction</b>                                      | <b>3</b>  |
| <b>Notices and disclaimers</b>                           | <b>4</b>  |
| <b>Trademarks</b>  | <b>4</b>  |
| <b>Principles of the test</b>                            | <b>5</b>  |
| <b>Kit contents</b>                                      | <b>8</b>  |
| <b>Reagents and equipment to be supplied by the user</b> | <b>8</b>  |
| <b>Kit storage</b>                                       | <b>9</b>  |
| <b>Suitable sample material</b>                          | <b>9</b>  |
| <b>Licensing agreement and limitations of use</b>        | <b>9</b>  |
| <b>Primerdesign satisfaction guarantee</b>               | <b>9</b>  |
| <b>Quality control</b>                                   | <b>9</b>  |
| <b>Bench-side protocol</b>                               | <b>10</b> |
| <b>qPCR amplification protocol</b>                       | <b>12</b> |
| <b>Interpretation of results</b>                         | <b>13</b> |

# Introduction

Janus kinase 2 (commonly called JAK2) is a human cytoplasmic tyrosine kinase protein that has been implicated in signalling by members of the type II cytokine receptor family (e.g. interferon receptors), the GM-CSF receptor family (IL-3R, IL-5R and GM-CSF-R), the gp130 receptor family (e.g. IL-6R), and the single chain receptors (e.g. Epo-R, Tpo-R, GH-R, PRL-R).

Mutations within the JAK2 gene are implicated in a wide range of myeloproliferative disorders and fusions with the TEL (ETV6) (TEL-JAK2) and PCM1 genes have been found in leukaemia patients. This mutation (V617F) results in a change of valine to phenylalanine at the 617 position and so renders hematopoietic cells more sensitive to growth factors such as erythropoietin and thrombopoietin.

Detection of JAK2 (V617F) is of diagnostic significance and quantification of this mutation is also useful in monitoring patients as a residual disease marker.

## Notices and disclaimers

During the warranty period Primerdesign quasa detection kits allow precise and reproducible data recovery combined with excellent sensitivity. For data obtained by violation to the general GLP guidelines and the manufacturer's recommendations the right to claim under guarantee is expired.

PCR is a proprietary technology covered by several US and foreign patents. These patents are owned by Roche Molecular Systems Inc. and have been sub-licensed by PE Corporation in certain fields. Depending on your specific application you may need a license from Roche or PE to practice PCR. Additional information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Roche Molecular Systems, 1145 Atlantic Avenue, Alameda, CA 94501 or Applied Biosystems business group of the Applied Biosystems Corporation, 850 Lincoln Centre Drive, Foster City, CA 94404. In addition, the 5' nuclease assay and other homogeneous amplification methods used in connection with the PCR process may be covered by U.S. Patents 5,210,015 and 5,487,972, owned by Roche Molecular Systems, Inc. and by U.S. Patent 5,538,848, owned by The Perkin-Elmer Corporation. The purchase of Biosearch Technologies products does not, either expressly or by implication, provide a license to use this or other patented technology. Licensing information can be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, CA 94404 or the Licensing Department at Roche Molecular Systems Inc., 1145 Atlantic Avenue, Alameda, CA 94501."

## Trademarks

Primerdesign™ is a trademark of Primerdesign Ltd.

Precision® is a registered trademark of Primerdesign Ltd.

genesig® is a registered trademark of Primerdesign Ltd.

The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

ABI, ABI PRISM®, GeneAmp® and MicroAmp® are registered trademarks of the Applied Biosystems Corporation.

BIOMEK® is a registered trademark of Beckman Instruments, Inc.; iCycler™ is a registered trademark of Bio-Rad Laboratories, Rotor-

Gene is a trademark of Corbett Research. LightCycler™ is a registered trademark of the Idaho Technology Inc.

GeneAmp®, TaqMan® and AmpliTaqGold® are registered trademarks of Roche Molecular Systems, Inc.,

The purchase of the Primerdesign™ reagents cannot be construed as an authorization or implicit license to practice PCR under any patents held by Hoffmann-La Roche Inc.

quasa A Master Mix containing GoTaq® Hot Start Polymerase manufactured by Promega Corporation for distribution by Primerdesign Ltd Licensed to Promega under U.S. Patent Nos. 5,338,671 and 5,587,287 and their corresponding foreign patents


# Principles of the test

## quasa (quantitative allele specific amplification)

quasa is a method for sensitive detection of mutations that may be present in clinical samples at low levels. quasa utilizes Primerdesign's own method based on 'allele specific PCR', using **modified primers** and hydrolysis probe, **modified cycling conditions** and a **modified master mix**. In this method the 3' terminal base of the mutation detecting primer is sited to bind to the mutant base. In the samples where the mutation is present, efficient amplification results in detection of the mutant sequence whereas PCR from this primer on Wild type samples is blocked due to the mismatch.


## Modified primers

quasa primers are designed such that the 3' terminal base overlies the mutation site. Thus, Wild type primers confer 100% specificity with the Wild type sequence but have a single base mismatch with the mutated sequences (The converse is true of the mutant specific primers). This is typical of allele specific PCR and relies on the principle that the single base mismatch will prevent the Wild type primers successfully priming on the mutated template. However, this principle alone is often insufficient to confer specificity and false positive amplification is common. Thus, quasa primers are further modified to multiply the level of specificity possible. Firstly, the primers are designed to have a low melting T<sub>m</sub> such that a single base mismatch will be as disruptive as possible. quasa primers also contain a sequence independent 'tag' at the 5' end. The tag is incorporated in to the amplicon during the first round of PCR and is thus present in the amplicon for subsequent cycles. This means that the tagged primers will prime preferentially on this template and thus drive amplification of the correct sequence in subsequent cycles. The quasa method requires no other primers, clamping primers, modified bases or blocking probes to achieve the extraordinary sensitivity of the kit.



## Modified cycling conditions

The quasa protocol uses a two-stage cycling strategy. The first five cycles of PCR use a very low annealing temperature of 50°C. This allows the low  $T_m$  primers to prime successfully whilst conferring the highest possible levels of specificity. After the first five cycles the annealing temperature is switched to 60°C. Thus, allele specific priming is effectively blocked and priming only occurs where the tagged primers have been incorporated. This also drives highly efficient amplification and probe cleavage and thus confers the sensitivity of the method. It is essential that the correct cycling parameters are used to ensure both sensitive and selective detection for each primer mix on its cognate target sequence.



## Modified master mix

Even greater levels of specificity and sensitivity are conferred by the specialized master mix used with quasa. Enzyme, Salt, dNTP and  $MgCl_2$  levels are all highly optimised. The annealing temperatures of the primer and probe have been carefully calibrated and any change in the reaction buffer can significantly alter the performance of the assay. For this reason, Primerdesign can only guarantee accurate results when Precision quasa A Master Mix is used.

## Quantitative analysis

The kit employs a two-reaction protocol. The user runs one reaction to detect the level of Wild type sequence present and at the same time another reaction to detect the level of mutant sequence present. By comparing the detection levels obtained by both primer/probe sets, the proportion of sequences that contain the mutation can be measured and expressed as a percentage relative to the Wild type sequences. When using the kit for a quantitative test each sample must be run in duplicate in order to get accurate results. The kit is sensitive down to detection levels of 0.1% of mutant DNA but at these levels the copy number of mutant DNA will always be very low and close to the limits of PCR detection. When working close to the limits of detection PCR becomes less accurate and it is therefore necessary to take the average of two points to ensure that a reliable result is calculated.

### **Qualitative uses of the kit**

In some clinical scenarios it is sufficient to know if a mutation is present or absent and the exact proportion of the mutant is not of diagnostic or therapeutic value. In these circumstances, the sample can be tested using a single test for each of the Wild type and mutant primer/probe mixes.

### **Positive template containing 1% mutant sequence**

The kit contains a positive control which contains a blend of both Wild type and mutant sequences at a known copy number. The ratio of Wild type to mutant template is 100:1 which is typical of some biological samples which can contain mutant sequences at a very low level. The positive control therefore provides a template for both primer and probe mixes and the quantification cycle (Cq) data from this control is used in the quantitative analysis.

### **Negative control**

To confirm absence of contamination, a negative control reaction should be included every time the kit is used. In this instance the RNase/DNase free water should be used instead of template. A negative result indicates that the reagents have not become contaminated while setting up the run. If a positive result is obtained the results should be ignored and the test samples repeated. Possible sources of contamination should first be explored and removed.

## Kit Contents

- **JAK2 V617F specific primer/probe mix (50 reactions BROWN)**  
FAM labelled
- **JAK2 WT specific primer/probe mix (50 reactions BROWN)**  
FAM labelled
- **Precision® quasa A Master Mix (2x 0.5ml ORANGE)**
- **Positive control template (JAK2 V617F 1% mutant) (RED)**
- **Template preparation buffer (YELLOW)**  
For resuspension of positive control template
- **RNase/DNase free water (WHITE)**

## Reagents and equipment to be supplied by the user

- **Real-Time PCR Instrument**  
Must be able to read fluorescence through the FAM channel
- **Pipettes and tips**
- **Vortex and centrifuge**
- **Thin walled 0.2 ml PCR reaction tubes**



## Kit storage

This kit is transported frozen and should be stored at -20°C on arrival. Freeze/thawing cycles should be kept to a minimum once resuspended. Under these conditions reagents are stable for six months from date of resuspension.

## Suitable sample material

All kinds of sample material suited for PCR amplification can be used. Please ensure the samples are suitable in terms of purity, concentration and DNA integrity. Always run at least one negative control with the samples. To prepare a negative control, replace the test sample with RNase/DNase free water.

## Licensing agreement and limitations of use

PCR is covered by several patents owned by Hoffman-Roche Inc and Hoffman-LaRoche, Ltd. Purchase of Primerdesign kits does not include or provide license with respect to any patents owned by Hoffman-La Roche or others.

## Primerdesign™ satisfaction guarantee

Primerdesign takes pride in the quality of all our products. Should this product fail to perform satisfactorily when used according to the protocols in this manual, Primerdesign will replace the item free of charge.

## Quality control

As part of our ISO9001 and ISO13485 quality assurance system, all Primerdesign products are monitored to ensure the highest levels of performance and reliability.

## Bench-side protocol

To minimise the risk of contamination with foreign DNA, we recommend that all pipetting be performed in a PCR clean environment. Ideally this would be a designated PCR cabinet. Filter tips are recommended for all pipetting steps. The positive control template is a significant contamination risk and should therefore be pipetted after negative control and sample wells.

### 1. Pulse-spin each tube in a centrifuge before opening.

This will ensure lyophilised primer and probe mix is in the base of the tube and is not spilt upon opening the tube.

### 2. Resuspend the primer/probe mix in the DNase/RNase free water supplied, according to the table below:

To ensure complete resuspension, vortex each tube thoroughly.

| Component – resuspend in water                            | Volume |
|---|--------|
| <b>Pre-PCR box</b>  |        |
| JAK2 V617F Genotyping primer/probe mixes ( <b>BROWN</b> ) | 60µl   |
| JAK2 WT Genotyping primer/probe mixes ( <b>BROWN</b> )    | 60µl   |

### 3. Resuspend the positive control template in the template preparation buffer supplied, according to the table below:

To ensure complete resuspension, vortex each tube thoroughly.

| Component – resuspend in template preparation buffer | Volume |
|--|--------|
| <b>Post-PCR heat sealed envelope</b>                 |        |
| Positive control template ( <b>RED</b> )*            | 500µl  |

\* This component contains high copy number template and is a VERY significant contamination risk. It must be opened and handled in a separate laboratory environment, away from the other components.

### 4. Prepare DNA detection mix for each of the SNP and Wild type primer mixes

Prepare sufficient reactions for all samples and extra reactions to include wells for the positive control and a 'no template control'.

| Component                                    | e.g. 1 reaction | e.g. 10 reactions |
|--|-----------------|-------------------|
| 2X Precision quasa A Master Mix              | 10µl            | 100µl             |
| Genotyping primer/probe mix ( <b>BROWN</b> ) | 1µl             | 10µl              |
| RNase/DNase free water (WHITE)               | 4µl             | 40µl              |
| <b>Final volume</b>                          | <b>15µl</b>     | <b>150µl</b>      |

**5. Dispense 15µl of each master mix according to your plate lay out**

**Sample plate layout for 5 samples for Qualitative analysis**

|                | Sample 1 | Sample 2 | Sample 3 | Sample 4 | Sample 5 | +ve ctrl | NTC |
|----------------|----------|----------|----------|----------|----------|----------|-----|
| JAK2 V617F     | ●        | ●        | ●        | ●        | ●        | ●        | ●   |
| JAK2 Wild Type | ●        | ●        | ●        | ●        | ●        | ●        | ●   |

**Sample plate layout for 5 samples for Quantitative analysis**

|                        | Sample 1 | Sample 2 | Sample 3 | Sample 4 | Sample 5 | +ve ctrl | NTC |
|------------------------|----------|----------|----------|----------|----------|----------|-----|
| JAK2 V617F             | ●        | ●        | ●        | ●        | ●        | ●        | ●   |
| <b>Duplicate wells</b> | ●        | ●        | ●        | ●        | ●        | ●        | ●   |
|                        |          |          |          |          |          |          |     |
| JAK2 Wild Type         | ●        | ●        | ●        | ●        | ●        | ●        | ●   |
| <b>Duplicate wells</b> | ●        | ●        | ●        | ●        | ●        | ●        | ●   |

**6. Dispense 5µl of extracted DNA sample according to your plate set up**

To obtain a strong signal whilst minimizing the addition of possible PCR inhibitors from the extraction, the recommended concentration of DNA is 5ng/µl (25ng in total). For highest detection sensitivity, more template can be added but the concentration should not exceed 10ng/µl (50ng in total). For the 'No Template Control' substitute 5µl of RNase/DNase free water.

NB. If you would like to add a larger volume of DNA (up to 9µl) simply reduce the volume of water at step 4. The final reaction volume is 20µl.

**7. Dispense positive control template.**

Add 5µl of diluted positive control template to the positive control wells. The positive control well will therefore automatically contain 100,000 copies of WT DNA and 1,000 copies of mutant DNA. The final reaction volume is 20µl.

# qPCR amplification protocol

## Recommended

The following cycling conditions must be used to ensure both sensitive and efficient detection of mutant sequences:

|            | <b>Step</b>            | <b>Time</b> | <b>Temp</b> |
|------------|------------------------|-------------|-------------|
|            | qPCR Enzyme Activation | 2 min       | 95°C        |
| x5 cycles  | Denaturation           | 10s         | 95°C        |
|            | Annealing              | 15s         | <b>50°C</b> |
|            | Extension              | 15s         | 72°C        |
| x40 cycles | Denaturation           | 10s         | 95°C        |
|            | Annealing *            | 30s         | <b>60°C</b> |
|            | Extension              | 15s         | 72°C        |

\*Fluorogenic data should be collected through the **FAM** channel.

# Interpretation of results

The percentage of SNP present in the sample is calculated using the delta Cq method. The proportions of SNP and WT in the sample are corrected by reference to a positive control standard where the SNP is present at a known proportion of 1%.

## Sample data and calculations

The calculation is performed in two stages. Firstly, the delta Cq values are used to calculate relative detection levels between the biological sample and the 1% control for both the Wild type and the mutant. These relative amounts are then converted into a percentage.

### Equations

JAK2 WT Relative Amount =  $[2^{\Delta} - (\text{JAK2 WT Sample} - \text{JAK2 WT Control})] \cdot \text{Control Proportion of WT}$


JAK2 V617F Relative Amount =  $[(2^{\Delta} - (\text{JAK2 V617F Sample} - \text{JAK2 V617F Control})) \cdot \text{Control Proportion of SNP}]$

JAK2 V617F Percentage conversion

$$= \frac{\text{JAK2 V617F Relative Amount}}{\text{JAK2 WT Relative Amount} + \text{JAK2 V617F Relative Amount}} \cdot 100$$


### 1. Calculate percentages for samples and control data

#### Example data:



#### Test sample

JAK2 Wild type Sample Cq = 26  
JAK2 Wild type Control Cq = 24



#### Positive control

JAK2 V617F Sample Cq = 32  
JAK2 V617F Control Cq = 32

Delta Cq calculation =  $2^{\Delta} - (26-24) = 0.25$   
Ratio calibration =  $0.25 \cdot 100 = 25$

Delta Cq calculation =  $2^{\Delta} - (32-32) = 1$   
Ratio calibration =  $1 \cdot 1 = 1$

## 1. Calculate the Percentages

From the example set of data the Wild type to mutant ratio is 25:1, when calibrated against the 1% positive control. Calculating the percentage of either the Wild type or the mutant can be calculated from this ratio. For the mutant, the percentage conversion is as follows

JAK2 V617F Percentage conversion

$$= \frac{1}{25 + 1} * 100$$

$$= \mathbf{3.85\%}$$

Interpretation: “3.85% of the sample DNA is mutated in background of 96.15% Wild type.”

**The kit is sensitive down to detection levels of 0.1%. Results that report detection at a lower proportion than 0.1% should be considered as negative.**